

Structural motif search across the protein-universe with Folddisco

Hyunbin Kim, Rachel Seongeun Kim, Milot Mirdita, and
Martin Steinegger

Structural Motifs

- **Short, recurring 3D arrangements** of tertiary structural elements
- Form **recognizable patterns** across diverse proteins
- Often linked to:
 - Stability
 - Binding interactions
 - Catalytic activity / active sites

Why Do They Matter?

- **Functionally constrained:** Evolution preserves them at sub-Ångström resolution
- Motif identification can reveal functional clues, even when:
 - No known homologs exist
 - The protein's function is uncharacterized


Motif-Function Examples

- **Cys₂-His₂ zinc finger**

Stabilizes **DNA-binding domains** in transcription factors

- **CWxP, NPxxY, DRY** motifs in **GPCRs**

Drive **receptor activation** and signal transduction



Gap in Annotation Tools

- Most methods rely on **sequence** → function relationship

*However, this relationship is **indirect***

- Sequence determines structure, but **function is executed by structure**
Similar sequences ≠ identical functions
(due to structural divergence, context, or local geometry)
- **Distant residues in sequence** can form **critical 3D motifs**
- This reliance stems from:
 - The wide availability of **high-throughput sequencing** and **alignment** tools
 - **Limited structural data** (until recently)
 - Historically **inefficient structure comparison** methods
- In contrast, methods that model **structure** → **function** relationships:
 - Can provide **more direct functional insights**
 - Especially useful for **motif-level functional prediction**


Structure-based functional annotation

Challenges with Conventional Tools

- **Scarcity of structural data** (compared to sequence data)
- **Slow and computationally intensive** structural alignment
- Poor scalability to large databases or motif-scale queries


Foldseek: A Scalable Alternative

- Converts 3D structure into a **1D sequence** using a **3Di alphabet**
- Enables **fast and scalable structural alignments**
 - Treats structure comparison as sequence alignment
- **Not designed for motifs**
 - Assumes **linear residue matching**, which is effective for **global or domain-level alignments**
 - **Structural motifs are non-linear**, they often involve:
 - **Non-contiguous residues** in sequence
 - **Spatially close but sequentially distant** fragments



RCSB motif search

- Addresses **non-linear motifs** by:
 - **Breaking proteins into proximal residue pairs**
- **Features for Each Pair:**
 - Residue 1 **AA identity**
 - Residue 2 **AA identity**
 - **C α –C α distance**
 - **C β –C β distance**
 - **Angle** between C α –C β vectors
- Stored as a **5-feature set** in an **inverted index**
- Index maps to: **PDB entry + positions**
- **Scales with residue count**
- Indexing requires **~75 \times more operations** than residue count
 - Due to pairwise feature extraction & storage



Limitations of motif search tools


- The indexing time and storage requirement
 - RCSB motif search took **3.5 days and 55GB** to index 160,467 structure
 - pyScoMotif took **20.5 hours** for 195,000 structures, but still required **73G**
 - a faster Python-based motif finder utilizing the same pair representation, except that it **uses side-chain centroids instead of C β atoms**
- **Lack of flexibility** in handling various query motif types and length
 - RCSB supports query motifs of up to **10 residues**
 - Alignment-based fragment search methods can **handle longer, discontinuous queries**, but **struggle with short motifs** like catalytic triads or zinc fingers

FoldDisco

- The first motif search algorithm that supports **both short motif queries and long, discontinuous segments** within a single framework.
- Massive scale efficiency: indexing **53M structures in under 24 hours** (<1.5 TB) with queries taking only a **few second**
- Folddisco examines proximal residue pairs
 - Extracts and encodes feature sets, storing them in an index
 - Builds upon **RCSB's feature set** with
 - **Torsion angles** (N–C α , C β) from **trRosetta**
 - Capture **side-chain orientation**

Folddisco is a fast tool for sensitive motif detection in millions of protein structures

- Given motif-defining query residues it examines **proximal pairs (<20Å)** and computes **feature sets for each pair**.
- Each set is encoded and rapidly searched against **a precomputed index of pairwise features** from database structures
- Extended search:** it can generate additional feature sets accounting for amino-acid substitutions, side-chain flexibility, and increased distances/angles.



FoldDisco: Feature Set & Indexing Strategy




Feature Set Construction

- Identify **proximal residue pairs**:
 - **5 RCSB features** (black):
 - **2 new features** (pink): Torsion angles (N–Ca and Cβ)
- Each **7-feature set** is **bit-encoded**
- Stored in an **index**:
 - Maps feature sets → structure **IDs** where they occur
 - **No need** to store residue positions

Motif Querying Process

- Apply **same feature extraction** on query motif's proximal residues
- Perform a **“pre-filter” step**:
 - Retrieve all structure **IDs** with matching feature sets
- Optionally:
 - Post-process retrieved structures Match their residues (pink) to the query (gray)

FoldDisco: Pairwise Features & Efficient Encoding



Encoding Feature Sets as 32-bit Integers

- **AA types** (20 options): 5 bits each (AA_1 , AA_2 : 10 bits)
- **Distances** (0–20Å, 16 bins): 4 bits each (Ca , $C\beta$: 8 bits)
- **Angles** (cos & sin, 4 bins each): 4 bits per angle ($Ca-C\beta$, torsion angles: 12 bits)
- Total: **30 bits** + 2 padding bits = **32-bit unsigned integer**
- Two feature sets per pair (both AA_1-AA_2 and AA_2-AA_1 directions since dihedral angles are **non-symmetric**).

Indexing Phase

- **Assign unique IDs** to each protein structure
- Identify **proximal residue pairs** (within 20Å radius)
- For each pair:
 - Extract **two sets of 7 features**
 - **Encode** each feature set as a **32-bit unsigned integer**
- Use each integer as a **key** in the index →
 - Maps to **structure IDs** where the feature set appear

FoldDisco: Querying

Querying Phase

- Extract proximal residue pairs from the **query motif**
- For each pair (i, j), compute:
 - Feature set for (i, j) and (j, i)
(due to **asymmetry** in dihedral features)
- Encode feature sets as 32-bit integers

Extended search (optional):

- More encodings for each query proximal pair in given range
- Allow **AA substitutions**
- Looser **distance/angle thresholds**

Pre-Filtering Step

- Use query integers as **keys** to retrieve matching **structure IDs**

Scoring & Ranking

- **Rank matches** by:
 - Number of shared feature sets
 - **Rarity** of those sets (higher rarity → higher

$$\text{IDF}_e = \log_2 \left(\frac{\# \text{ total structures}}{\# \text{ structures with encoding } e} \right)$$

- **Coverage score**
 - Measures how well a candidate structure **covers** the query motif
 - Adjusted by structure length:

n: number of shared encodings

L: structure length (residues)

α: length penalty exponent (default = 0.5)

$$\text{Score}_{\text{cand}} = L^{-\alpha} \sum_{i=1}^n \text{IDF}_{e_i}$$

- **Motif Completeness Score**
 - Counts **distinct query residues** involved in shared encodings (e.g. (x-y and x-z) or (x-y and z-t))

FoldDisco: Residue matching via Graph Construction

Why Needed?

- FoldDisco's index doesn't store residue positions: Must match residues post hoc for **structural alignment**

Residue Matching Process

- Build **residue graph**:
 - Nodes = **candidate residues**
 - Directed edge if a residue pair matches any **query pair encoding**
 - Edges may also be added for similar **AA identity** and **Ca–Ca distance**
 - **More than two feature sets** are considered for each query residue pair by setting **the distance and angle thresholds**
- **Graph Search**
 - Identify **motif-like residue clusters** as:
 - **Strongly Connected Components** (via **Tarjan's Algorithm**)
 - **Weakly Connected Components** (via DFS on undirected graph)
- **Superposition computation**
 - Superposes the query motif on the matched residues using the Quaternion Characteristic Polynomial algorithm.
 - RMSD is calculated using the coordinates of the Ca and C β atoms of the query motif and the matched residues.

Folddisco is the most accurate method in querying the human fraction of the AFDB-proteome for zinc fingers, both when using a short motif query suitable for pyScoMotif and RCSB (d, left; residue labels, e.g. F207, denote chain and residue number) and when using the motif-containing segments suitable for MASTER (e, left). f, Folddisco achieves higher sensitivity than pyScoMotif on SCOPe-constructed benchmarks, where the goal is to match SCOPe sequences of the same family as the query before matching a different fold, using all conserved columns ("full") or a random subsample of them (60%, 20%).

Folddisco accurately detects discontinuous motifs like zinc fingers and segment-based motifs, previously requiring separate tools.

Dataset

Human subset of AFDB: 23,391 protein structures





Compared Against

- **Short motif:**
 - RCSB Motif Search
 - pyScoMotif
- **Segment-based motif**
 - MASTER

Zinc Finger Motif PDB: 1G2F

Full motif: F207, F212, F225, F229

Segment for MASTER: F204–215 and F222–232



Evaluating FoldDisco's Generalizability with SCOPe Benchmark

Performance beyond specific motifs (e.g., zinc fingers, catalytic triads)

Benchmark Construction


- Based on **SCOPe family-level MSAs** (from **FoldMason**)
- Selected **fully occupied columns** and a **dominant residue** (occurring in >66% of the members)

Simulates **realistic, scattered motif-like queries**

Three Query Types using the dominant residues:

- **Full**: All dominant residue positions
- **60% Sampled**: Random subset of positions
- **20% Sampled**: Sparse queries with minimal info

Folddisco achieves higher sensitivity than **pyScoMotif**




Evaluation Criteria

- **TP**: Match from same **SCOPe family**
- **FP**: Match from a **different fold**
- **Sensitivity**: TP / P before first FP

Ranked by:

- **Coverage score** (FoldDisco pre-filter)
- **RMSD** (FoldDisco full)


Folddisco builds indexes faster and smaller than previous tools:



indexing AFDB50 (53M structures) takes only ~24h vs. ~20 days (extrapolated) for pyScoMotif. Querying a zinc-finger motif across AFDB50 takes just ~13s, up to 48x faster than pyScoMotif.

Applications of Folddisco: Zinc finger motif detection

Folddisco can annotate proteins: querying a canonical zinc-finger uncovers an **uncharacterized oyster protein and metagenomic proteins**. It also detects **partial catalytic metal sites in E. coli** peptide deformylase. All of these hits would be **missed by Foldseek or sequence aligners**




Applications of Folddisco: Conformational state identification

Folddisco can distinguish functional states.

- Searching GPCR activation motifs (CWxP, NPxxY, DRY), clearly separating active/inactive states.
 - from activated (left, magenta)
 - inactivated (right, purple)

Large-Scale Search on AFDB


- ~**53%** of retrieved structures were **active**
- Closely aligns with **experimental PDB distribution** (~**54%** active)



Applications of Folddisco: Protein interface search

Folddisco queried a cross chain protein–protein interface motif pattern derived from immunoglobulin λ -like and immunoglobulin κ variable domains

- an interface between antibody chains (gray/black), it successfully identifies matching interfaces within monomeric antibody fragments (cyan).





Folddisco Server


- Databases
 - PDB
 - AFDB
 - ESM30


⌵


Folddisco Search



 Search


 Multimer search


 FoldMason MSA

 Folddisco search

History

< >

 2025-07-21 12:52
jv_3rgxbcb1PCF1xa...

 2025-07-18 19:10
XFD7QSKmFTF1obJ3Pd...

Results for job: XFD7QSKmFTF1obJ3PdHdhrPuLyU0rMaan3W0Vtg

ALL DATABASES

AFDB-PROTEOME (1000)

AFDB50 (1000)

ESM30 (1000)

PDB100 (1000)

AFDB-PROTEOME 1000 hits

SHOW TAXONOMY

Filter

Query residues

F207

F212

F225

















F229

Cluster

☐ Cluster

Clustering Min Points
2

Clustering Epsilon
8

Target	Description	Scientific Name	idf-score	RMSD	Nodes	Matched residues	Structure
AF-Q86VK4-F1-MODEL_V4	Zinc finger protein 410	Homo sapiens	5.511	0.167	4	A251,A256,A269,A273	
AF-A0A0R4IX80-F1-MODEL_V4	Sik1-like domain 2	Danio rerio	8.669	0.176	4	A20,A25,A38,A42	
AF-P22227-F1-MODEL_V4	Zinc finger protein 42	Mus musculus	7.761	0.184	4	A229,A234,A247,A251	
AF-Q8ST83-F1-MODEL_V4	Polcomb protein PHO	Drosophila melanogaster	6.755	0.186	4	A416,A421,A434,A438	
AF-B4F792-F1-MODEL_V4	Zinc finger protein 410	Rattus norvegicus	6.259	0.187	4	A251,A256,A269,A273	
AF-D4A6Z4-F1-MODEL_V4	Zic family member 4	Rattus norvegicus	7.448	0.193	4	A198,A203,A216,A220	
AF-E7F8R1-F1-MODEL_V4	Zinc finger protein 410	Danio rerio	7.197	0.193	4	A221,A226,A239,A243	
AF-A0A1S805M9-F1-MODEL_V4	Uncharacterized protein	Dracunculus medinensis	5.706	0.198	4	A380,A385,A398,A402	
AF-O1L884-F1-MODEL_V4	WT1 transcription factor b	Danio rerio	4.997	0.201	4	A318,A323,A336,A340	
AF-F1Q751-F1-MODEL_V4	Uncharacterized protein	Danio rerio	4.127	0.202	4	A220,A225,A238,A242	
AF-Q54ET8-F1-MODEL_V4	C2H2-type zinc finger-containing protein	Dictyostelium discoideum	4.016	0.203	4	A598,A603,A616,A620	
AF-A0A175WEP1-F1-MODEL_V4	Uncharacterized protein	Madurella mycetomatis	9.382	0.203	4	A18,A23,A36,A40	
AF-Q8BMU0-F1-MODEL_V4	Zinc finger protein 76	Mus musculus	8.862	0.203	4	A257,A262,A275,A279	
AF-F1M349-F1-MODEL_V4	Zic family member 5	Rattus norvegicus	8.419	0.206	4	A95,A100,A113,A117	
AF-A0A5K4F2Z1-F1-MODEL_V4	Uncharacterized protein	Schistosoma mansoni	6.080	0.207	4	A619,A624,A637,A641	
AF-A0A0K9G075-F1-MODEL_V4	Uncharacterized protein	Stenodrilus stenosilis	5.609	0.209	4	A539,A544,A557,A561	

Conclusion

High-Speed Motif Search

- Indexes **millions of structures** in under 24 hours
- Queries return results in **seconds**, with high sensitivity

Motif-Centric and Functionally Aware

- Handles both **short/local** and **long/discontinuous** motifs
- Distinguishes **functional states** (e.g., active vs. inactive GPCRs)
- Enables **interface-level searches** for **PPI discovery**

Beyond Global Alignment

- Moves past sequence and domain-based matching
- Captures **structural motifs** across diverse folds and conformations